

# STALWART

**Eyepiece Adapter (Reduction Lens)  
SCN-0.45X-2**



## Description

These adapters are used to connect the C-mount cameras to the microscope eyepiece tube or trinocular tube of 23.2mm. If the eyepiece tube diameter is 30mm or 30.5mm, you can plug the 23.2 adapter into the 30mm or 30.5mm connecting ring and then plug into the eyepiece tube.

## Specification

	<b>C-mount</b>	<b>Magnification</b>	<b>Connecting Dia.</b>	<b>Application</b>
SCN0.45x-1 adapter	Yes	0.45x	23.2mm	Used to connect camera to eyepiece / trinocular tube
SCN0.45x-2 adapter	Yes	0.45x	23.2mm	Used to connect camera to eyepiece / trinocular tube
SCN0.5x adapter	Yes	0.5x	23.2mm	Used to connect camera to eyepiece / trinocular tube
SCN30 Connecting Ring	No	No	30mm	Used with adapters to plug into 30mm eyepiece tube
SCN30.5 Connecting Ring	No	No	30.5mm	Used with adapters to plug into 30.5mm eyepiece tube