

STALWART

DSLR Camera Adapter
SDCA-Series



Description

These 2 adapters are used to connect DSLR camera to the microscope eyepiece tube or trinocular tube of 23.2mm. If the eyepiece tube diameter is 30mm or 30.5mm, you can plug the 23.2 adapter into the 30mm or 30.5mm connecting ring and then plug into the eyepiece tube.

Specification

	Camera Mount	Magnification	Connecting Dia.	Application
SDCA-1 (NIKON) adapter	Nikon	2x	23.2mm	Used to connect Nikon DSLR camera to eyepiece / trinocular tube
SDCA-2 (CANON) adapter	Canon	2x	23.2mm	Used to connect Canon DSLR camera to eyepiece / trinocular tube