

STALWART

Laboratory Water Bath STW-22 Series



Introduction

Suitable for direct heating and auxiliary heating of biological, chemical, physical, plant, chemical and other experimental samples. Routine laboratory temperature control, Escherichia coli detection, sample thawing, bacterial detection, incubation microbial experiments, cell culture, food testing pretreatment, etc.

Features

- PID temperature controller, large LCD display screen and interface provides for user- friendly operation.
- Easy to set Adjustable Timer. (1 minute to 5,999 minutes)
- Built-in circulating water pump to ensure uniform upper and lower temperature of bath lotion. (only for BW-22P)
- Maintenance-free operation with easy to clean surface.
- Both heater and bath chamber are made of corrosion-resistant stainless steel.
- No angle in bath chamber, easy to clean.
- The standard stainless-steel bottom plate, helps prevents direct contact by accessories and tubes to heating element and sensors.
- Drain valve makes emptying of bath water fast and easy for cleaning and moving.

Safety

- Preset On/Off function.
- Built-in power interruption protection function, automatic run after power interruption.
- Independent Over-temperature protection meets DIN 12880 International standard requirements.
- Temperature deviation alarm.
- Over current protection alarm.

Option

- Available RS-485 or USB ports for data collection.

Specification

Item	STW-22	STW-22P
Power requirements	AC220V 50HZ	
Power (W)	1000W	
Temperature range	RT+5 ~ 100°C	
Temperature Fluctuation	±0.2°C	
Temperature Resolution	0.1°C	
Chamber volume	22L	
Internal Dimension W×D×H(mm)	280×220×150	
External Dimension W×D×H(mm)	345×285×355	
Timer	1 ~ 5999min	
Net Weight	18 Kg	19 Kg
Porous cover(lid)	6-hole	6-hole

Note: STW-22P is with built-in circulating water pump